i · ϕ · s[®]
flexible fixation

IOS
Instructions
For Use

| | |
|------------------|----|
| English | 1 |
| Deutsch | 2 |
| Français | 3 |
| Español | 4 |
| Italiano | 5 |
| Nederlands | 6 |
| Portugués | 7 |
| Ελληνικά | 8 |
| Svenska | 9 |
| Norsk | 10 |
| Dansk | 11 |
| Suomi | 12 |
| 日本語 | 13 |

Instructions for use

Intended use

INTELLIGENT ORTHOPAEDICS instruments consist of manual surgical instruments and devices intended for use in surgical procedures.

This instrument is a **Class II single use device for use with Intelligent Orthopaedics Fracture Management systems ONLY.**

Instrument utilisation is determined by the user's experience and training in surgical procedures. Do not use this instrument for any purpose outside the intended use of the device, as it may seriously affect the safety and function of the product.

CAUTION: Handle devices with care to prevent damage to surgical gloves.

Recommendations for sterilisation

INTELLIGENT ORTHOPAEDICS recommends that before sterilisation the products are decontaminated and cleaned. The cleaning and decontamination of instruments and devices must follow the guidelines set forth by MHRA, AORN/HIMA and AAMI. Both physical and chemical (detergent) processes are necessary to minimise the bioburden on all soiled items. Chemical (detergent) cleaners alone cannot remove all soil and debris therefore a careful manual cleaning of each item with soft sponge or cloth is essential for maximum decontamination. Carefully inspect hidden areas, such as recesses, to assure any residual materials are removed. Once the items have been cleaned and decontaminated they should be thoroughly rinsed with clean water to remove any detergent or chemical residue before sterilisation. The use of mild enzymatic detergent is recommended.

Materials used in INTELLIGENT ORTHOPAEDICS devices can be sterilised using steam sterilisation methods. MHRA, AORN/HIMA, and AAMI guidelines for sizes, weights and mass should be followed. This product has been certified sterilisable using the following cycle: Porous Load 134-137°C / 3 minutes.

1. After sterilisation all instruments should be allowed to cool. The time of cooling is dependent on the load size and mass. Place instruments on a rack or shelf with linen cover until cooling is complete. The potential for condensation may increase if the cases are not allowed to cool properly.

CAUTION: Hot instruments should not be handled or used as they can cause injury.

2. If condensation is observed check to ensure that cooling, as described in 1. above, has been followed. Verify that the steam being used for the sterilisation process has a quality more than 97%. Also confirm that the sterilisers have been inspected for routine maintenance in accordance with manufacturer's recommendations.

3. Suggested steam exposure times: (ref: HTM2010 Part 3)

Preferred: 134-137°C / 3 minutes
Alternatives: 126-129°C / 10 minutes
121-124°C / 15 minutes
115-118°C / 30 minutes

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©2007 BSN medical Inc.
Intelligent Orthopaedics, the IO logo, the IOS logo and the storm logo are all trademarks of Intelligent Orthopaedics Ltd.
BSN medical Ltd.
Intelligent Orthopaedics Ltd.